



since 1990



UNIQUE AND POWERFULL

CONTINUOUS INHALATION

DRY SALT THERAPY



INDIVIDUAL CONTINUOUS INHALATION WITH(OUT) MASK, ALSO INTENDED FOR HORSE TRAILER

ULTRASONIC NEBULIZER

Technical Data

- Nebulization rate: 0.0-3.0 ml/min
- Air volume: max 20 l/min
- Particle size: 0.5-5.0 microns
- Capacity of medication cup: max 150 ml (optimal 50 ml)
- Ultrasonic frequency: 1.7 MHz
- Noise level: <35 dB (A)
- Operating voltage/frequency: 100-240 V AC, 50/60 Hz or 12V DC
- Power consumption: 55 W
- External dimensions: 280 x 200 x 220 mm
- Weight of the main unit: 3.3 kg



Contents of the kit:

PROFI SONIC A ultrasonic nebulizer, inhalation hose (90 cm), 3 medication cups (one already installed), medication supply set, AC power cord, DC power cord and instruction manual.

Additional option: Equine mask

ULTRASONIC ATOMIZER

Technical Data

- 3 piezoelectric crystals, which are glass coated due to resistance to chemicals
- Capacity of atomization: 0-1.2 l/h
- The most common particle size: 5.0 microns
- Ultrasonic frequency: 1.7 MHz
- Power supply/frequency: 100-240 V AC, 50/60 Hz
- Power consumption: 150 W
- Dimensions: 320 x 280 x 270 mm
- Weight: 7,2 kg



PRIZnano HAS A CAPACITY
OF TRANSFORMING 1200ml
OF LIQUID INTO A LARGE
AMOUNT
OF AEROSOL
WITHIN ONE
HOUR

CE





GROUP OR INDIVIDUAL DRY SALT THERAPY,
WITHOUT DRUGS, FOR HORSE TRAILERS OR
LARGE INDOOR SPACES

DRY SALT AEROSOL GENERATOR

Technical Data

- Touch screen operated
- Personalization of up to 9 programs:
duration of inhalation and salt dosing
- Capacity of the salt container: max 40 g
- Salt consumption: 0-100 g/h (adjustable)
- Operating voltage/frequency:
100-240 V AC, 50/60 Hz
- Power consumption: 18 W
- Dimensions: 225 x 222 x 403 mm
- Weight: 8 kg

multiSALT+ IS A COMPACT
DEVICE WHICH GENERATES
DRY AEROSOL FINE, DRUG
FREE PARTICLES, FOR
EFFICIENT HALOTHERAPY





since 1990

PRIZMA IS LIMITED COMPANY WITH A **TRADITION OF 31 YEARS**

WE PRODUCE, DISTRIBUTE AND SELL **MEDICAL EQUIPMENT**

WE HAVE BEEN THE PRODUCER SINCE 1998

OUR PRODUCTION AND SALES RESPECT STANDARDS:

ISO13485 / ISO 9001 / ISO 14001

OUR PRODUCTS QUALITY IS CERTIFIED BY **CE / UL MARKS**





Prizma GmbH

GERMANY

40213 Düsseldorf, Graf-Adolf-Platz, 15

Tel. +49 211 97632861

e-mail: s.katanic@prizma.gmbh

SERBIA

34000 Kragujevac, N° 8, Kumanovska st.

Tel: +381 34 330 200

Fax: +381 34 336 146

e-mail: prizma@prizma.rs

www.prizma.rs