



since 1990



UNIQUE AND POWERFULL

CONTINUOUS INHALATION

DRY SALT THERAPY



## INDIVIDUAL CONTINUOUS INHALATION WITH (OUT) MASK, ALSO INTENDED FOR HORSE TRAILER

### ULTRASONIC ATOMIZER

#### Technical Data

- Atomization rate: 0.0-3.0 ml/min
- Air volume: max 20 l/min
- Particle size: 0.5-5.0 microns
- Capacity of medication cup: max 150 ml (optimal 50 ml)
- Ultrasonic frequency: 1.7 MHz
- Noise level: <35 dB (A)
- Operating voltage/frequency: 100-240 V AC, 50/60 Hz or 12V DC
- Power consumption: 55 W
- External dimensions: 280 x 200 x 220 mm
- Weight of the main unit: 3.3 kg



#### Contents of the kit:

PROFI SONIC A ultrasonic atomizer, inhalation hose (90 cm), equine mask, 3 medication cups (one already installed), medication supply set, AC power cord, DC power cord and instruction manual.





## ENVIROMENTAL INHALATION FOR HORSES - NO NEED FOR ADDITIONAL ACCESSORIES

### ULTRASONIC ATOMIZER

#### Technical Data

- 3 piezoelectric crystals, which are glass coated due to resistance to chemicals
- Capacity of atomization: 0-1.2 l/h
- The most common particle size: 0.5-5.0  $\mu\text{m}$
- Ultrasonic frequency: 1.7 MHz
- Power supply/frequency: 100-240 V AC, 50/60 Hz
- Power consumption: 150 W
- Dimensions: 250 x 250 x 350 mm
- Weight: 6,5 kg

PRIZnano HAS A CAPACITY  
OF TRANSFORMING 1200ml  
OF LIQUID INTO A LARGE  
AMOUNT  
OF AEROSOL  
WITHIN ONE  
HOUR





GROUP OR INDIVIDUAL DRY SALT THERAPY,  
WITHOUT DRUGS, FOR HORSE TRAILERS OR  
INDOOR SPACES

## DRY SALT AEROSOL GENERATOR

### Technical Data

- Touch screen operated
- Personalization of up to 9 programs:  
duration of inhalation and salt dosing
- Capacity of the salt container: max 35 g
- Salt consumption: 0-100 g/h (adjustable)
- Operating voltage/frequency:  
100-240 V AC, 50/60 Hz
- Power consumption: 18 W
- Dimensions: 225 × 395 × 218 mm
- Weight: 8 kg

**multiSALT+ IS A COMPACT  
DEVICE WHICH GENERATES  
DRY AEROSOL FINE, DRUG  
FREE PARTICLES, FOR  
EFFICIENT HALOTHERAPY**





since 1990

PRIZMA IS LIMITED COMPANY WITH A **TRADITION OF 33 YEARS**

WE PRODUCE, DISTRIBUTE AND SELL **MEDICAL EQUIPMENT**

WE HAVE BEEN THE PRODUCER SINCE 1998

OUR PRODUCTION AND SALES RESPECT STANDARDS:

**ISO13485 / ISO 9001 / ISO 14001**

OUR PRODUCTS QUALITY IS CERTIFIED BY **CE / UL MARKS**





Prizma GmbH  
GERMANY  
40213 Düsseldorf, Graf-Adolf-Platz, 15  
Tel. +49 211 97632861  
e-mail: s.katanic@prizma.gmbh  
[www.prizmade.de](http://www.prizmade.de)

SERBIA  
34000 Kragujevac, N° 8, Kumanovska st.  
Tel: +381 34 330 200  
Fax: +381 34 336 146  
e-mail: [sales@prizma.rs](mailto:sales@prizma.rs)  
[www.prizma.rs](http://www.prizma.rs)