




Prizma GmbH
GERMANY
40213 Düsseldorf, Graf-Adolf-Platz, 15
Tel. +49 211 97632861
e-mail: s.katanic@prizma.gmbh
www.prizmade.de

Prizma Ltd
SERBIA
34000 Kragujevac, N° 8, Kumanovska st.
Tel: +381 34 330 200
Fax: +381 34 336 146, 336 147
e-mail: sales@prizma.rs
www.prizma.rs/en



Health Care

Professional



since 1990



Our Story

SINCE 1990

1990
year that we were founded

Production Experience
26 years

2019, foundation of Prizma GmbH,
based in Germany

EU, USA, Indonesia, Australia
our export markets

From 1990 Omron Healthcare Distributor
Serbia, Bosnia & Herzegovina, Montenegro,
and North Macedonia markets

Cooperation with
Rwth Aachen University, Serbian Academy of Science

CE  CE 1304
our products certificates

ISO 9001:2015, EN ISO 13485:2016, ISO 14001:2015
our management system standards

MEDICAL DEVICES

breathe care



PRIZMA

EFFICIENT OXYGEN DELIVERY SYSTEM

PRIZMA O2, portable oxygen concentrator, provides a high oxygen concentration to patients, with oxygen needs.

The device is delivering medical-grade oxygen through the nasal cannula, while removing nitrogen.

When a patient needs oxygen, a system for effective and efficient oxygen delivery is essential.

PRIZMA O2 is a lightweight device that comes with a battery and can be used at home, at work or while travelling.

PRIZMA O2 delivers oxygen through a pulse dose using intelligent technology.

PRIZMA O2 has simple control functions and an LCD screen.



DEVELOPED & MADE IN SERBIA

PRIZMA O2

CE 1304*

- **Oxygen delivery on the go**
- **Lightweight and quiet**
- **Easy to use and easy to maintain**



Technical Data

- Oxygen Concentration: 87-96%
- Flow type: pulse
- Levels of flow settings: 6
- Flow rate: 210-1260 ml/min
- Max limited pressure: 1.5 bar \pm 25%
- Battery life: up to 5 h, at 210 ml/min
- Battery charging time: 2.5 h
- Noise level: \leq 60 dB
- Power Supply: 100-240 V AC, 50-60 Hz; 24 V DC, 120 W
- Warm up time: 2 min
- Dimensions: 187 \times 88 \times 217 mm
- Weight: 2.16 kg

CONTENTS OF THE KIT:

PRIZMA O2 oxygen concentrator, 8-core battery, AC/DC adapter, nasal cannula, carrying bag and instruction manual

* Q1 2025

PRI/MA

PROFESSIONAL ULTRASONIC NEBULIZERS

Prizma is producing the own line of professional ultrasonic nebulizers for medical facilities; devices are simple to use and easy to maintain.

Our professional ultrasonic nebulizers are ideal for continuous multipatient nebulization.

On the device's touch screen it is possible to make settings such as: air volume, nebulization rate, start time and duration of inhalation therapy. Screen is readable even in the dark, thanks to the backlight.

Our product range consists of:

PROFISONIC – basic professional ultrasonic nebulizer

PROFISONIC H – professional ultrasonic nebulizer with integrated aerosol heating system



DEVELOPED & MADE IN SERBIA

PROFI SONIC

CE 1304

- Large number of patients continuous inhalations
- Touch screen for parameters setting
- In every position nebulization



Technical Data

- Nebulization rate: 0.0-3.0 ml/min
- Air volume: max 20 l/min
- Particle size: 0.5-5.0 microns
- Capacity of medication cup: 5-150 ml
- Low drug residue: 1 ml
- Capacity of water tank: 300 ml
- 100% silicone hose:
ultra-thin silicone hose can be sterilized in an autoclave
- Ultrasonic frequency: 1.7 MHz
- Noise level: <35 dB (A)
- Operating voltage/frequency: 100-240 V AC, 50/60 Hz
- Power consumption: 55 W
- Dimensions: 260 × 200 × 250 mm
- Weight: 1.8 kg

CONTENTS OF THE KIT:

PROFI SONIC nebulizer, silicone inhalation hose (90 cm), mouthpiece, 3 medication cups (one already installed), adult mask, child mask, bacterial filter, power cord and instruction manual.

PROFI SONIC H

CE 1304

- **Heated nebulization**
continuous inhalations
- **Adjusting the heat level**
of nebulization
- **Touch screen**
for parameters setting

Technical Data

- Nebulization rate: 0.0-3.0 ml/min
- Air volume: max 20 l/min
- Particle size: 0.5-5.0 microns
- Capacity of medication cup: 5-150 ml
- Low drug residue: 1 ml
- Capacity of water tank: 300 ml
- 100% silicone hose:
ultra-thin silicone hose can be sterilized in an autoclave
- Ultrasonic frequency: 1.7 MHz
- Noise level: <35 dB (A)
- Operating voltage/frequency: 100-240 V AC, 50/60 Hz
- Power consumption: 90 VA
- Dimensions: 260 × 200 × 250 mm
- Weight: 1.8 kg

CONTENTS OF THE KIT:

PROFI SONIC H nebulizer, silicone inhalation hose (90 cm), 3 medication cups (one already installed), bacterial filter, power cord and instruction manual.

Optional accessories are available.



OPTIONAL ACCESSORIES

PROFI SONIC / PROFI SONIC H

Table stand, adjustable height
Made of stainless steel
Catalogue code: PSDS



PROFI SONIC / PROFI SONIC H

Movable stand, adjustable height
Made of stainless steel
Catalogue code: PSTS



PROFI SONIC / PROFI SONIC H

Medication supply set
Catalogue Code: PSMS
Material: PP (polypropylene)
Single pack quantity: 1 piece



PROFI SONIC / PROFI SONIC H

Inhalation hose
Catalogue Code: PSIH
Material: Silicone
Single pack quantity: 1 piece



PROFI SONIC / PROFI SONIC H

Bacterial filter
Catalogue Code: PSBF
Material: Polypropylene
Single pack quantity: 1 piece



PROFI SONIC / PROFI SONIC H

Medication cup
Catalogue Code: PSMC
Material: PET
Single pack quantity: 5 pieces



PROFI SONIC

Mouthpiece
Catalogue Code: PSMP
Material: PC
Single pack quantity: 5 pieces



PROFI SONIC

Child mask
Catalogue Code: PSCM
Material: Silicone
Single pack quantity: 1 piece



PROFI SONIC

Adult mask
Catalogue Code: PSAM
Material: Silicone
Single pack quantity: 1 piece



PRI/MA

COMPRESSOR NEBULIZER

Prizma is producing powerful compressor nebulizer with unique nebulizer kit that allows setting of nebulization rate of 0.3-0.45 ml/min, providing efficient drug delivery to the smallest parts of the airways and making the therapeutic effect of the administered drug even greater.

During nebulization the users are protected of inhaling dust and other impurities from the air, thanks to the filter that is supplied in our compressor nebulizer.

Our compressor nebulizer is intended for use in the home environment and medical facilities, and patients could be adults and kids.



DEVELOPED & MADE IN SERBIA

PRIZJet

CE 1304



- **Setting of aerosol rate with unique nebulizer kit**
- **Strong compressor**
- **For home & medical facilities use**

Technical Data

- Particle size: <5 microns
- Nebulization rate: 0.3-0.45 ml/min
- Capacity of medication cup: 2-8 ml
- Low drug residue: 0.64 ml
- Compressor maximum pressure: 2.5 bar
- Compressor maximum air flow: 8 l/min
- Noise level: <60 dB (A)
- Operating voltage/frequency: 230 V/50 Hz
- Power consumption: 120 VA
- Dimensions: 180 × 140 × 110 mm
- Weight: 1.3 kg (compressor only)

CONTENTS OF THE KIT:

PRIZJet nebulizer, the power cord, inhalation hose (length 150 cm), nebulizer kit (medication cup), adult mask, child mask, mouthpiece, air filters (5 pcs + 1 assembled), storage and carrying case pouch and instruction manual.

Optional accessories available for different types of inhalation therapy.



high precision

MEDICAL DEVICE



PRI/MA

ANEROID SPHYGMOMANOMETER

Prizma is producing aneroid sphygmomanometer with EC certificate.

100% German mechanism made by HB-Feinmechanik GmbH & Co. KG is the heart of our aneroid sphygmomanometer, that provides accurate measurements and reliable readings.

Its one-tube system guarantees easiest operation while unique designed 180°rotary anti-slip bulb enables both left- and right-hand users to easily use the device.

Bladderless cuffs are biocompatible according to the high criteria of EN ISO 10993 series of standards and produced with high quality TPU coated nylon fabric that are latex-free, not toxic, non phthalate.

It is intended for use in the medical facilities and in home environment.



DEVELOPED & MADE IN SERBIA

PRI/MA PA1

CE 1304



- **100% German mechanism** by HB-Feinmechanik
- **180° rotary anti-slip bulb** for left-handed and right-handed users
- **Bladderless cuff** biocompatible (EN ISO 10993)

Technical Data

- Manometer scale range: 0-300 mmHg
- Accuracy: ± 3 mmHg
- Scale diameter: $\varnothing 66$ mm
- Tube system: one-tube system
- Cuff type: bladderless for arm circumference 27.5-36.5 cm
- Cuff material: TPU coated nylon
- Measuring mechanism, German-made by HB-Feinmechanik GmbH & Co. KG


CONTENTS OF THE KIT:

PRI/MA PA1 aneroid sphygmomanometer, bladderless cuff 27.5-36.5 cm, storage and carrying case pouch and instruction manual. Optional cuffs for arm circumferences of 9.0-14.8 cm, 13.8-21.5 cm, 20.5-28.5 cm, 27.5-36.5 cm and 35.5-46.0 cm are available on demand



long tradition

INDUSTRIAL DEVICE



PRIZMA

R&D SCALE ULTRASONIC SPRAY PYROLYSIS SYSTEM

PRIZnano is a large capacity ultrasonic atomizer of liquids. Usage of 3 piezo-electric crystals guarantee a large atomization capacity up to 1.2 l/h, which provides production of large quantity of aerosol that is ideal for independent use (humidification of large spaces) or in combination with other system for the production of nano particles.

PRIZMA ESP nano particle collector use high-voltage electrostatic field to collect ultra-fine particles. The minimum collectable particle is 0.1 nm. Collection efficiency is more than 95%.

It is intended for use in industry, science and research & development centers (in the process of obtaining nano powders), in medicine and food industry.



DEVELOPED & MADE IN SERBIA

PRIZnano

CE



- **Pyrolysis atomizer**
our specialization
- **Resistant to chemicals**
the chamber
- **RWTH Aachen university**
developed with

Technical Data

- 3 piezoelectric crystals, which are glass coated due to resistance to chemicals
- Capacity of atomization: 0-1.2 l/h
- The most common particle size:
 - 5.0 microns – 1.7 MHz
 - 3.0 microns – 2.4 MHz
- Ultrasonic frequency: 1.7 or 2.4 MHz
- Power supply / frequency:
100-240 V AC, 50/60 Hz
- Power consumption: 150 W
- Dimensions: 250 × 250 × 350 mm
- Weight: 6.5 kg

CONTENTS OF THE KIT:

PRIZnano atomizer and instruction manual.

PRIZMA ESP

CE



- Easy collection of nano powder

- Easy cleaning



Collected silver nano particles

Technical Data

- Power supply / frequency:
100-240 V AC, 50/60 Hz
- Power consumption: 500 W
- Dimensions: 600 × 550 × 560
- Weight: 32.5 kg


CONTENTS OF THE KIT:

PRIZMA ESP filter and instruction manual.



Prizma GmbH
GERMANY
40213 Düsseldorf, Graf-Adolf-Platz, 15
Tel. +49 211 97632861
e-mail: s.katanic@prizma.gmbh
www.prizmade.de

Prizma Ltd
SERBIA
34000 Kragujevac, N° 8, Kumanovska st.
Tel: +381 34 330 200
Fax: +381 34 336 146, 336 147
e-mail: sales@prizma.rs
www.prizma.rs/en



Halotherapy

Dry Salt Aerosol Generators



since 1990



Halotherapy

What is dry salt therapy

What is Halotherapy?

Today, more and more people are suffering from respiratory issues. From pollution and smog to airborne diseases, pollens and toxins; never before has our society been faced with such a growing epidemic of respiratory ailments.

Halotherapy (dry salt therapy) has been researched over the past few decades.

WHAT IS DRY SALT THERAPY (HALOTHERAPY)?

Halotherapy (dry salt therapy) is a natural, drug-free, alternative treatment that may help to improve respiratory issues and skin conditions.

The key process involves a halogenerator crushing and grinding 99.99% Pure Grade Sodium Chloride (salt) into very tiny microparticles. The halogenerator, or halotherapy machine, then disperses these dry aerosol microparticles of salt into the salt therapy room, salt cave, halotherapy sauna, or salt booth during a typical session.

As an individual undergoes a salt therapy session, they inhale the microparticles. The particles penetrate deep into the lungs, bronchi, bronchioles, and alveoli to address respiratory concerns.

At the same time, the non-inhaled salt particles are absorbed by the body to combat skin conditions. Salt therapy can offer both respiratory (lungs and breathing) and cutaneous (skin-related) benefits. Additionally, the relaxing, calm environment offered by salt room therapy can help with stress reduction.

Most people are familiar with “wet” salt therapy. Whether you’ve put salt in a bath, had a saline solution treatment, used a Nettie pot to clear your nasal passages, or gargled with saltwater, wet salt has its benefits.

Dry salt room therapy is different. The drier the salt, the better.

That is why halotherapy protocols include heating the salt prior to putting it into the salt chamber. Heat ensures that all moisture is evaporated. When the salt is crushed and ground into microparticles and inhaled, the salt particles can do their best job.

The Benefits Of Salt Room Therapy

Halotherapy benefits adults and children alike, as well as athletes and animals. It is natural and safe and there are no side effects.

A majority of people who visit salt rooms on a regular basis may find relief from a variety of respiratory conditions. These include, among many others:

- Asthma
- Allergies
- Common cold
- Bronchitis
- COPD
- Cystic fibrosis
- Sleep Apnea

Dry salt therapy is great for people who also have particular skin conditions such as eczema, acne and psoriasis. Salt therapy may also provide help when someone is recovering from cosmetic surgery.

However, one does not need to have any specific conditions to gain benefits from dry salt therapy.

Halotherapy is very beneficial for overall general wellness. It removes toxicity from the respiratory system. People are exposed on a daily basis to pollutants, airborne diseases, bacteria, allergens and other irritating and harmful factors.

Halotherapy is an excellent way to cleanse and detox the lungs. It can also help to invigorate the whole body by supporting increased lung capacity and oxygen intake.

Interested in salt therapy for general respiratory hygiene? Visiting your local salt therapy facility once or twice a week will enhance lung capacity, boost the immune system, reduce stress, and support better sleep.

Many chronic and acute respiratory conditions require increased frequency and consistency of sessions to realize the most benefits. We encourage you to make salt therapy a habit and a regular part of your health and wellness routine.



Salt Room Therapy

The Benefits of Salt Room

PRIZMA

DRY SALT AEROSOL GENERATORS

No halotherapy without a halogenerator. Dry salt therapy (halotherapy) can only be provided by utilizing a professional-quality halogenerator. Salt particle size emitted from the halogenerator must be 5 microns or less (the smaller the better) Dry salt therapy (Halotherapy) can only be provided in an enclosed space to contain the salt particles in the area such as a room, booth, cabin, or other type of chamber.

Prizma is producing the own line of dry salt aerosol generators with EC and UL certificates for medical facilities and wellness / SPA centers.

On the device's touch screen or via the Android / Windows / IOS smart device is possible to make settings of device, such as salt consumption, start time, duration of halotherapy, dosage profile (device is offering 10 dosage profiles). Screen is readable even in the dark, thanks to the backlight.

Materials, which are in contact with salt, are made of stainless steel and plastic, that ensures long life of Prizma dry salt aerosol generators.

Our product range consists of:

- multiSALT+** – wellness / SPA device with UL certificate and CE mark for small salt rooms, salt cabins and tent
- PRIZSALT+** – medical device with EC certificate
- PRIZSALT+S** – wellness / SPA device with UL certificate and CE mark
- homeSALT** – home device with UL certificate and CE mark to be used also in salt room, tent, cabin or in outdoor spaces
- nanoSALT** – home device with UL certificate and CE mark to be used also in salt room, tent, cabin or in outdoor spaces



DEVELOPED & MADE IN SERBIA

multiSALT+



- Wellness/SPA device
- Wireless management Prizma Link app
- Up to 9 programs personalization



Technical Data

- Touch screen operated
- Personalization of up to 9 programs: duration of inhalation and salt dosing
- Capacity of the salt container: max 35 g
- Salt consumption: 0-100 g/h (adjustable)
- Operating voltage/frequency: 100-240 V AC, 50/60 Hz
- Power consumption: 18 W
- Dimensions: 225 × 218 × 395 mm
- Weight: 8 kg

CONTENTS OF THE KIT:

multiSALT+ dry salt aerosol generator, cleaning equipment, power cord, deflector, extension tube, software and instruction manual.



PRIZSALT+

CE 1304



- **Medical device**
- **Wireless management**
Prizma Link app
- **Up to 9 programs**
personalization

Technical Data

- Touch screen operated
- Personalization of up to 9 programs:
duration of inhalation and Salt dosing
- Capacity of salt dispenser: max 150 g
- Salt consumption: 0-500 g/h (adjustable)
- Operating voltage/frequency:
100-240 V AC, 50/60 Hz
- Power Consumption: 35 W
- Dimensions: 270 x 260 x 480 mm
- Weight: 13.5 kg

CONTENTS OF THE KIT:

PRIZSALT+ dry salt aerosol generator, cleaning equipment, power cord, deflector, extension tube, software and instruction manual.



PRIZSALT+S

UL CE



- **Wellness/SPA device**
- **Wireless management**
Prizma Link app
- **Up to 9 programs**
personalization

Technical Data

- Touch screen operated
- Personalization of up to 9 programs:
duration of inhalation and Salt dosing
- Capacity of the salt dispenser: max 150 g
- Salt consumption: 0-500 g/h (adjustable)
- Operating voltage/frequency:
100-240 V AC, 50/60 Hz
- Power consumption: 35 W
- Dimensions: 270 x 260 x 480 mm
- Weight: 13.5 kg

CONTENTS OF THE KIT:

PRIZSALT+S dry salt aerosol generator, cleaning equipment, power cord, deflector, extension tube, software and instruction manual.



homeSALT



- **New approach in halotherapy**
- **Home device for indoor and outdoor spaces**
- **Remote control**
4 programs of operation

Technical Data

- Remote control: Operate from a distance
- Programs of operation: 4 programs of operation
- Led indicator: Light indication of device status
- Operating voltage: 24 V
- AC/DC adapter
100-240 V AC, 47-63 Hz, 1.6-0.7 A
24 V DC, 2.6 A max
- Power consumption: 15 W
- El.protection class: Class III
- Salt tank capacity: 20 g maximum
- Dimensions: 113 × 236 × 70 mm
- Net weight: 1 kg

CONTENTS OF THE KIT:

homeSALT dry salt aerosol generator, AC/DC adapter, remote control key fob, brush and instruction manual.

nanoSALT



- **New approach in halotherapy**
- **Small dimension battery operated**
- **Portable**

Technical Data

- 10 inhalations with a fully charged battery
- Battery charge indicator
- Empty battery indicator
- Indicator of blocked motor
- Open cover indicator
- Operating voltage: 6V DC
- AC/DC adapter:
100-240 V AC, 50-60 Hz, 0.5 A
5 V 3 A DC
- Power consumption: 8 W
- El.protection class: Class III
- Salt tank capacity: 5 g
- Dimensions: 170 × 50 × 60 mm
- Net weight: 0.5 kg

CONTENTS OF THE KIT:

nanoSALT dry salt aerosol generator, AC/DC adapter, brush and instruction manual.



The following authorities spoke about the benefits of halotherapy

during STAcon LIVE 2024 in Orlando Florida, in occasion of 10th anniversary of Salt Therapy Association

Allan Share President of the Spa Industry Association



As President of the Spa Industry Association with thousands of connections across the globe and a dedicated and knowledgeable team, I have learned some valuable lessons about what it takes to succeed in the spa and wellness world. Now more than ever is the time to be sharing these insights with one another, and so, I'll continue doing what I enjoy and share some of this information.

Leo Tonkin Founding Director of the STA



Leo has been instrumental in developing salt/halotherapy into an entire wellness modality and industry. He initiated establishing the Salt Therapy Association 10 years ago working with the leading researchers, physicians, healthcare providers, manufacturers, and owners. Twelve years ago, he started SALT Chamber, the pioneering company leading the salt therapy and wellness industry with over 3,600 projects and installations. He brings over 30 years of leadership, culture, and strategy consulting expertise to entrepreneurial start-ups to FORTUNE 100 companies working with Stephen Covey, Ken Blanchard and others. He leverages his extensive experience by collaborating with other industry experts in his role as the founding chairman of the Salt Therapy Association and Chair of the Global Wellness Institute's Respiratory Wellness Initiative. He actively shares his knowledge by presenting at conferences like Biohacking, Cryotherapy, Floatation, and others, engaging with various sectors of the global wellness industry, including fitness, hospitality, and spas. He is also a cancer survivor.

Dr. Mark Erath Professor Mayo Clinic. TED-TALK Presenter



Dr. Mark Erath, M.D.: Physician, Scientist, Inventor, Social Entrepreneur, and TEDx Speaker. Emeritus professor at Mayo Clinic College of Medicine & Science. During his 25-year tenure at Mayo, Dr. Erath authored over 100 peer-reviewed publications and presented more than 250 lectures globally. For the past decade, he has focused on combating airborne disease transmission and addressing the health impacts of air pollution. In this field, he has evaluated various technologies, worked as a consultant, and delivered a TEDx Talk titled "One Less..." on April 6, 2022, at TEDx University.

Jean Guy de Gabriac Founder of World Wellness Weekend



Award-winning educator in Guest experience, specializing in elevating skills and engagement for Hotel and Spa teams. Celebrating 20 years as CEO of TIP TOUCH International (2004-2024). Founder of World Wellness Weekend, celebrated by 8,000+ venues in 150 countries with 60 supporting organizations. The 8th edition is scheduled for 20-22 Sept 2024 (9th edition: 19-21 Sept 2025). Additional roles include Conference producer for World Spa & Wellness Convention London, Board Member of the International Massage Association, and judge for both the World Massage Championship and World Spa & Wellness Awards.

Shoshana Treichel Above and Beyond Yoga & Salt Therapy Center



After a 1986 attack in NYC, Shoshana became a world-record-holding powerlifter and top women's bodybuilder. She owns Above and Beyond Yoga and Salt Therapy facility, develops university weightlifting curricula, and trains elite athletes. As a trainer, nutritionist, and yoga instructor, she's helped clients transform their bodies, working with NFL, NBA, and MMA professionals.

Debra Hope Duneier EchoChi



An esteemed Master Feng Shui practitioner, accredited LEED® Green Associate, Certified Eco-Designer, and Graduate Gemologist, Debra will be leading a session on "Designing Wellness Experiences & Environments". Learn how to create spaces that enhance well-being and sustainability.

Steve Woods Salt Mind Health / Promerica Health



Former NBA and MLB agent, tennis promoter, and sports marketing professor. Co-owner of an NBA team and Owner of Restore Hyper Wellness group. CEO of Promerica Health, a global health engagement company. He is also the founder and owner of Salt Mind Health.

Craig Ure Marketing & Business Growth



Digital marketing strategist with almost three decades of experience in business growth strategy and implementation. Turning this experience and specialized skillset into digital transformations and automations for SMB's leveraging AI solutions and architecture. Host of the STA bi-weekly Marketing Mastermind.

Peter Kater

Multi-Grammy winning and Multi-Platinum selling pianist, composer, producer and breathwork practitioner



Kater will lead an immersive exploration of how music and breathwork can transform well-being. Attendees will experience the power of breathwork firsthand in a guided session from his new SPIRITUS breathwork album, followed by a live Q&A session, providing a unique opportunity to engage with the artist. This session promises to be a profound experience, highlighting the therapeutic potential of combining music and breathwork for holistic health. The term "Spiritus" is Latin for "breath" or "to breathe," reflecting the profound connection between breathwork and holistic wellness.

Dr. Tami Peavy MBA, MPT, DPT, CEO CEO and Clinical Director of La Mesa Rehab



Dr. Tami Peavy is a distinguished Doctor of Physical Therapy with over 20 years of experience and advanced degrees from the Medical College of Georgia, Northwestern University, Indiana University, and Utica University. As the CEO and Clinical Director of La Mesa Rehab in San Diego, she specializes in pulmonary rehabilitation, offering a unique program that combines respiratory therapy and physical therapy treatments. Dr. Peavy leads a multidisciplinary team providing high-quality care for patients.

Eddie Fisher

Chairman of The Board for The American Lung Association of Orlando & Southeast Regional Board



Financial advisor with 23+ years experience in personalized wealth management strategies. Chairman of the Orlando Florida Board and Southeast Regional Board of Directors for the American Lung Association. Received 2022 Southeast Region Board Member of the Year. Supports lung health causes in memory of his parents.

Dr. Daniel Layish Board Certified Pulmonary Disease, Critical-care



Dr. Layish is a pulmonologist with experience in a wide range of lung conditions, including COPD, pleural effusion, and sleep apnea. While he treats many respiratory issues, his particular areas of interest include pulmonary hypertension, sleep disorders, cystic fibrosis, alpha-1 antitrypsin deficiency, and lymphangioleiomyomatosis.

Tara Grodjesk

President of Earthlite Spa and Wellness Products



As a consultant to resorts and hotels in the creation of new and exciting signature spas including products and treatment program menus. Tara has trained over 3000 massage therapists and aestheticians internationally and continues to educate spa staffs in the wisdom and techniques of Aromatherapy and Ayurveda.

Emmanuel Verduzco Co-Founder / CEO - Ageless Fitness Corp.



CEO and co-founder of Ageless Fitness Corp. Serial entrepreneur with 4 startups focused on services for older adults. Former social worker in California. Launched a homecare franchise in 2014, followed by ventures in non-emergency medical transport, specialized fitness for seniors, and advocacy. Successfully franchised his adaptive fitness program for older adults in 2020.

Amy Davis

Owner of Just Breathe Salt Therapy Wellness Spa



Amy, owner of Just Breathe Salt Therapy Wellness Spa in Winter Park, FL, transitioned from a 32-year TV sales career to wellness. She applies her marketing expertise to grow the spa, offering innovative salt therapy treatments. Amy enjoys travel and beach time with family, including her two children and five grandchildren.

Clayton Joyner Founder and CEO of MODRN Digital



Clayton Joyner, CEO of MODRN Digital, is a leader in digital growth strategies. With executive experience as CGO and CMO, he guides businesses through digital transformation. His expertise in marketing technology helps companies navigate the complex digital landscape, driving innovation and optimizing growth in the modern era.

Brayden D. Whitlock, JD, PhD
 Director of Research Outbreaker Solutions



Former biotech professional turned lawyer with a PhD in toxicology. Specializes in facilitating technical collaborations and transactions. Experienced educator skilled at explaining complex concepts. Serves on various boards and committees. Teaches health innovation, law, and research methodology. Published researcher and advocate before courts, tribunals, and regulators including the US EPA and Health Canada.

Jimmy Youngblood
 President - Osteostrom /
 Inventor of the Red Halo FX



A pioneering figure in biohacking for nearly a decade, Jimmy serves as the President and CTO of Osteostrom. He is a visionary inventor at the forefront of cutting-edge health and wellness technologies, with a focus on skeletal strength and red light therapy products. Renowned for harnessing the frontiers of science to drive unprecedented breakthroughs, Youngblood is the driving force behind technologies that redefine bone health and vitality. His innovations hold immense financial promise for business owners and entrepreneurs. Jimmy has been a pioneer in developing red light therapy and is the inventor of the "Red Halo FX," recognized as the most versatile and powerful system in its field.

Lars Lofstrand
 Cryotherapy Industry National Association



A "Serial Entrepreneur" who has started, owned, and run over 22 different businesses. Passionate about helping aspiring and established entrepreneurs achieve success. A respected leadership authority, author, teacher, speaker, futurist, and organizational consultant. Bestselling author of "Sticky Money" and "Real Estate Game Changers." Consistently rated as a top speaker and facilitator at global conferences and seminars. His appeal transcends age, culture, and occupational barriers.

Tammy Pahel
 VP of Spa & Wellness at Carillon Miami



Tammy's tenure at Carillon Miami has been marked by a commitment to elevate spa and wellness operations. They've successfully integrated technology and medically focused treatments to set new industry benchmarks. Empowering guests to achieve optimal health, and their team is driven by the collective goal of nurturing individual potential.



Be part of our network, join us!

Distribution Network

- **Europe**
 - **Slovenia** CST PLUS d.o.o.
Krivce 5 1000 Ljubljana
 - **Italy** BIORESEARCH SRL
Palazzo Manzoni - Via dei Mille 14 25122 Brescia
 - **Romania** SC HIMALAYASALT SRL
Str. Intr. Victor Daimaca nr.6, sector 1, 13279 jud. Bucuresti
 - **Hungary** FOCUS HOUSE KFT
Taliga street 29 1031 Budapest
 - **Slovakia / Czech Republic** DASOLA spol. s.r.o.
Bakošova 22 841 03 Bratislava
 - **Montenegro** SPA MEDICA DOO
8. marta 74 81000 Podgorica
 - **UK / Iceland** SERENE SALT ROOMS
262 High Road, Harrow, England, HA3 7BB
- **North America**
 - **USA • Canada** S.A.L.T. Chamber LLC.
FL., 120 NW 11th Street, 33432 Boca Raton, USA